

How to Track Progress

HappyNeuron Pro tracks each exercise completed by a patient and provides visualization tools to help you follow-up.

Benefits of tracking

The HappyNeuron Pro program automatically tracks all your patients' results, whether they are with you or completing workouts at home. So it is really useful to see how they engage in the program. You can also use this information to show your patients and their caregivers how they are progressing. It can also help you to see where they are struggling and adapt your therapy accordingly.

Visualize the patient's results

To visualize a patient's results, you will need to complete the following steps:

Log into your account > Navigate to the "Users" section at the top of the screen > Under the user list, click on the patient whose results you would like to visualize. > Click on the "Results" tab




Patient Results Overview

This screen shows you all of the activity that was completed by your patient during the month. It allows you to see their scheduled workouts, the daily activity and their overall strength and weakness.

To change the month to visualize with the "Previous" and "Next" buttons.


You can easily print a printed version by selecting this button.




Breaking Apart the Exercises

To view more a in-depth of a particular exercise:


Click any exercise in the bottom overview chart



Click on "Overview" on the top right of the overview screen to open a drop down menu and choose an exercise.





In this Exercise View, choose a date range on the top and you will be able to visualize all the results on this exercise. You will see the number of games played and the current level.



Use the filter button to show only the results that matches your settings filter.

The first graph shows the level evolution of the patient in this exercise. Roll over or click any dot to get more in-depth details.



The second graph shows the progression of accuracy and time for this exercise.

Roll over or click any item to get details.

The last table shows a list of the results for this exercise, ordered by date.

	Accuracy	Speed	Variant Code	
Thursday, March 09, 2023				
	50%	5s	1 - 21111111	▶ +
	50%	7s	1 - 21111111	▶ +
	100%	8s	22121111	▶ +
Monday, March 13, 2023				
	50%	3s	1 - 21111111	▶ +
	100%	8s	4 - 22221111	▶ +
	50%	7s	1 - 21111111	▶ +

Here you can replay exercises with the same variables by clicking on this variable ▶ or you can add them to be completed in the workout tab by clicking on this variable + .

Patient Access

By default, patients who log in to their account have access to the **same results visualization module**, so they can also see their progression and engagement.

You can decide not to show this module to certain patients if you think it can be counter-productive. Just uncheck the option ["View Performance Profile"](#) in the patient file.

View Performance Profile

Allow user to choose variants

Access to exercises outside of workouts: